AH/AHS-DRY Range Automatic autoclaves with drying

AH-21 L

AHS-50 DRY

AH/AHS-DRY Range

S Class: Suitable to sterilize: unwrapped and wrapped solids, small porous objects and hollow ware type B (slightly deep and broad cavities), glass, plastic, metallic objects, waste bags.

AH-21 L is suitable to sterilize liquids.

Operating Principal

AH/AHS-DRY table-top range provides a single solution for the multiples sterilization needs on laboratories including instruments, hollowware elements, glassware and other laboratory items that need a final drying to avoid contamination.

Load has to be placed into the vessel's baskets and, after automatic filling with distilled water, device starts to heat up and purge until set temperature is reached. In case of solid objects, discharge is automatically forced during cooling stage to reduce cycle time.

Digital Microprocessor

Digital microprocessor for an easy programming and parameters selection for sterilization process.

Programs – Microprocessor has 10 programs, from P0 to P9, and the first four are pre-defined and protected.

Program N°	Temperature (°C)	Pressure (bar)	Sterilization time (min.)	Sterilization mode	Drying time (min.)
P0	115	0,6	60	SOLIDOS-1	12
P1	121	1,1	30	SOLIDOS-1	25
P2	133	2,1	20	SOLIDOS-1	30
P3	121	1,1	20	LIQUIDOS-1	-

The rest of the programs are user-defined with the following parameters setting:

- Sterilization temperature
- Sterilization time
- Drying time
- Heart probe selection

AH/AHS-DRY Range

Water Filling – Pure water for sterilization can be added manually into the water pan. Automatic filling is also available by connecting the device to pure mains water.

Sterilization Process

Pre-vacuum and filling stage – Integrated vacuum pump for initial pre-vacuum and subsequent automatic water filling.

Heating Stage – Fast saturated steam production by means of powerful heating elements.

Sterilization Control – A PT-100 temperature probe controls the chamber temperature, with the possibility of installing another flexible PT-100 temperature probe for liquid programs.

Cooling and drying stage – In solid program, discharge is automatically forced during cooling stage to reduce the cycle time. Then, integrated heating jacket together with pre-vacuum pulses is activated for a fast and intense drying.

Security

All safety standards are covered with safety valve, safety thermostat and safety pressure switch, door opening prevention system in case of overpressure, open door sensor and water level detector. EU Directive and Standards:

2014/35/UE *Low Voltage*

2014/30/UE *Electromagnetic Compatibility*

2014/68/UE *Pressure Equipment*

EN-61010-1

EN-61010-2-040

EN-61326

AD 2000 Merkblatt

Communications and Documentation


RS-232 Port – For live monitoring of sterilization cycle directly to PC. It is needed the customized software SW7000 (optional).

Integrated Printer (optional) – For quality printouts after each cycle. Intervals can be adjustable. Printer has to be requested before ordering.

AH/AHS-DRY Range


Components supplied with the autoclave AH-21 L

1. Tray support for four trays.
2. Three stainless steel perforated trays. Ref. BAH-21
3. Auxiliary tray.
4. A holding clamp.
5. 1.5 m flexible tube to drain condensation reservoir.
6. Bag support. Ref. BAP-21

Components supplied with the autoclaves AHS-50/75 DRY

1. Support for five trays.
2. Two perforated trays made of stainless steel. Ref. BAH-B
3. Auxiliary tray.
4. 1.5 m flexible tube to drain condensation reservoir.

AH/AHS-DRY Range

Technical Data

Referencia Reference Référence	Volumen total / útil Total volume / usable Volume total / utile	Ø x H mm	L x D x H mm	Temp. °C	Presión máx. bar Max. pressure bar Pression max. bar	W	Kg	Rango tiempo min. Time range min. Rang de temps min.	Depósito H ₂ O l Container H ₂ O l Réservoir H ₂ O l
AH-21-L	22/21	250 x 430	560 x 660 x 425	100-134	2,1	2000	55	1' - 250'	6
AHS-50 DRY	55/50	400 x 400	805 x 805 x 650	100-134	2,1	2800	125	1' - 250'	15
AHS-75 DRY	79/75	400 x 600	805 x 1005 x 650	100-134	2,1	3200	140	1' - 250'	20

AHS-50/75 models with closing wheel

Electrical connection for equipment of 21, 50 and 75 L: 230V (1P + N), 50/60 Hz.
Possibility of modifying the voltage on request.

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FLEXIBLE "HEART" PROBE PT-100

- Recommended for sterilizing liquids in flasks of 1 litre or bigger.
- The control is done from the microprocessor.
- This accessory should be installed in factory.

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THERMAL PAPER PRINTER

- It prints program number, cycle number, temperature, date and hour of the run.
- Error messages are also printed.
- To be installed in our facilities.

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Referencia: **ITS**

TABLE TOP PRINTER

- It prints program number, cycle number, temperature, date and hour of the run.
- Error messages are also printed.
- It includes RS232 connection.

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