

Digital thermostatic baths



BAD-4

Features:

Temperature regulation by microprocessor with digital read out and P.I.D. action. PT100 probe.
 Programmable delayed start (0-99,5 hours). Programmable incubation time (0-99,5 hours).
 Safety under Regulation EN 61010-1
 Safety thermostat with manual rearms according to DIN 12877 Class 2
 Inner tank made in AISI 304 stainless steel.
 External case made in AISI 304 Stainless steel heating elements inside the tank.
 Incorporated drainage

Series	Regulation	Temperature range	Temperature stability	Homogeneity	Resolution
BAD	Digital	Room temp. +5° - 99,9° C	± 0,15° C	± 1° C	0,1° C
BOD	Digital	Room temp. +5° - 200° C	± 0,4° C	± 2° C	1° C

Reference	Capacity l.	Int. H x W x D in mm	Ext. H x W x D in mm	Power W	Weight kg.
BAD-2	5	120 x 295 x 140	210 x 345 x 190	1000	5
BAD-4	12	120 x 295 x 320	210 x 345 x 470	1500	7,5
BAD-6	20	120 x 490 x 290	210 x 345 x 645	2000	10
BAD-12	40	120 x 620 x 500	210 x 665 x 645	3000	15
BOD-2	5	120 x 295 x 140	210 x 345 x 190	1000	5
BOD-4	12	120 x 295 x 320	210 x 345 x 470	1500	7,5
BOD-6	20	120 x 490 x 290	210 x 345 x 645	2000	10
BOD-12	40	120 x 620 x 500	210 x 665 x 645	3000	15

Digital thermostatic baths Accessories

Concentric rings lid

TB-6



Roof type lid

TT-6



Heater protecting tray

TI-6



Support frames

MB001, MB002, MB003



Bath series	Capacity l.	Heater protecting tray	Concentric rings lid	Roof type lid	Support frame
BAD, BOD	5	TI-2	TB-2	TT-2	MB001
BAD, BOD	12	TI-4	TB-4	TT-4	MB002
BAD, BOD	20	TI-6	TB-6	TT-6	MB003
BAD, BOD	40	TI-12	TB-12	TT-12	MB003 (2 u)

Rack with handles



Reference	Description
GRB-5613	56 tubes (14 x 4) Ø 13 mm
GRB-4816	48 tubes (12 x 4) Ø 16 mm
GRB-3020	30 tubes (10 x 3) Ø 20 mm
GRB-5011	50 tubes (10 x 5) Ø 11 mm